

Portable Kahuna Shoulder Brace



The Portable Kahuna Shoulder Brace is a universal shoulder brace that can be applied to either the right or left side of the patient.

HOW TO USE

1. Remove the Kahuna Shoulder Brace from the bag and separate the **Sling** from the **Abduction bar**.
2. On the **Abduction bar**:
 - a. Unclip the Velcro strap at the front of the Abduction bar.
 - b. Bend the larger, padded, malleable Abduction bar around the torso and waist.
 - c. Connect Velcro strap to the front of the Abduction bar and tighten like a belt.



- d. The secondary, shorter abduction bar, on the outside of the main abduction bar, will remain unbent and point forward (it should be perpendicular to the torso).



3. On the **Sling**:

- a. Unbuckle the Sling Strap from the front of the sling.
- b. Unstrap the two Velcro straps and slide arm into the sling.
- c. Grab onto the stress ball at the front-end of the sling.
- d. Secure arm inside the sling by re-strapping the two Velcro straps.
- e. Run the Sling Strap diagonally up your back, over your opposite shoulder, and then along the front of your torso.
- f. Connect the Sling Strap buckle to the front-end of the sling.



4. With the Abduction bar secured to the torso, and the arm secured inside of the Sling: Secure the Sling to the Velcro portion of the shorter abduction bar, which is on the outside of the main Abduction bar.

